

Human THO complex subunit 4 (THOC4) ELISA Kit

Catalog No: #EK20582



Package Size: #EK20582-1 48T #EK20582-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Human THO complex subunit 4 (THOC4) ELISA Kit
Brief Description	ELISA kits
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.312-20 ng/mL
Sensitivity:0.125ng/mL
Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm
Detection Method:sandwich

Note: This product is for in vitro research use only