

# Mouse S100 Calcium Binding Protein A8 (S100A8) ELISA Kit



Catalog No: #EK20306

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Package Size: #EK20306-1 48T #EK20306-2 96T

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Mouse S100 Calcium Binding Protein A8 (S100A8) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

## Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:35pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only