

Human S100 Calcium Binding Protein A8 (S100A8) ELISA Kit



Catalog No: #EK20305

Orders: order@signalwayantibody.com

Package Size: #EK20305-1 48T #EK20305-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human S100 Calcium Binding Protein A8 (S100A8) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:26pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Note: This product is for in vitro research use only