

USP32 Polyclonal Antibody

Catalog No: #49189

Package Size: #49189-1 50ul #49189-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	USP32 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IHC;IF;WB;ELISA
Species Reactivity	Human;Mouse
Specificity	USP32 Polyclonal Antibody detects endogenous levels of USP32 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human USP32.
Target Name	USP32
Other Names	USP32;USP10;Ubiquitin carboxyl-terminal hydrolase 32;Deubiquitinating enzyme 32;Renal carcinoma antigen NY-REN-60;Ubiquitin thioesterase 32;Ubiquitin-specific-processing protease 32
Uniprot	Q8NFA0
GeneID	84669
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

Images



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using USP32 Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only