

Human OGDHL/ 2-oxoglutarate dehydrogenase-like, mitochondrial ELISA Kit

Catalog No: #EK21134

Package Size: #EK21134-1 96T

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Human OGDHL/ 2-oxoglutarate dehydrogenase-like, mitochondrial ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	2-oxoglutarate dehydrogenase complex component E1-like, Alpha-ketoglutarate dehydrogenase-like, OGDC-E1-like, KIAA1290
Uniprot	Q9ULD0

Application Details

Detect Range:78-5000pg/mL

Sensitivity:32.5pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only