

# PerCP-Cy5.5 anti-mouse CD185

Catalog No: #26844



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PerCP-Cy5.5 anti-mouse CD185
Brief Description	Rat monoclonal
Host Species	Rat
Clonality	Monoclonal
Isotype	Rat IgG2b, k
Applications	FC
Species Reactivity	Mouse
Storage	Store at 4°C in dark

## Application Details

FC 5 $\mu$ /time

Note: This product is for in vitro research use only