

## NKG2D Rabbit Polyclonal Antibody Cy3 Conjugated

Catalog No: #C92910Cy3



Package Size: #C92910Cy3-1 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	NKG2D Rabbit Polyclonal Antibody Cy3 Conjugated
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	FC
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human NKG2D: 101-200/219
Conjugates	Cy3
Target Name	KLRK1
Other Names	KLRK1; natural killer cell group 2D; NKG2-D isoform b; Klrk1; Nkrp2; NKG2-D type II integral membrane protein; NKG2-D-activating NK receptor; NK cell receptor D; Killer cell lectin-like receptor subfamily K; NK lectin-like receptor; CD314; CD 314; CD314 antigen; D12S2489E; Killer cell lectin like receptor subfamily K member 1; Killer cell lectin-like receptor subfamily K member 1; KLR; Klrk1; NKG2 D activating NK receptor; NKG2 D type II integral membrane protein; NKG2-D-activating NK receptor; NKG2D_HUMAN; NKLLR; NKR P2; Nkrp2.
Uniprot	P26718
Concentration	1 mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

## Application Details

Flow-Cyt 1ug/Test

Note: This product is for in vitro research use only