

# HIF-1-alpha Antibody

Catalog No: #57174

Package Size: #57174-1 50ul #57174-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	HIF-1-alpha Antibody
Host Species	Rabbit
Clonality	Monoclonal
Isotype	Rabbit IgG
Purification	Affinity-chromatography
Applications	WB,ICC/IF
Species Reactivity	Human,Mouse,Rat
Immunogen Description	A synthesized peptide derived from human HIF-1-alpha
Other Names	HIF1; MOP1; PASD8; bHLHe78; HIF-1alpha; HIF1-ALPHA; HIF1A
Calculated MW	110kDa
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at +4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.


## Application Details

WB 1:500~1:2000 ICC/IF 1:50~1:200

## Images



All lanes use the Antibody at 1:1K dilution for 1 hour at room temperature.



Western blot analysis of HIF-1-alpha expression in HeLa treated with CoCl<sub>2</sub> cell lysate.

---

Note: This product is for in vitro research use only