

VPS52 rabbit pAb

Catalog No: #49241



Package Size: #49241-1 50ul #49241-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	VPS52 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	This antibody detects endogenous levels of VPS52 at Human, Mouse,Rat
Immunogen Description	Synthesized peptide derived from human VPS52
Target Name	VPS52 SACM2L
Uniprot	Q8N1B4
GeneID	6293
Calculated MW	80kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only