

Canine Insulin (INS) ELISA Kit

Catalog No: #EK20712



Package Size: #EK20712-1 48T #EK20712-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Canine Insulin (INS) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Canine

Application Details

Detect Range: 3.12-200ng/mL

Sensitivity: 1.16ng/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time: 1-4.5h

Detection wavelength: 450 nm

Product Description

Detection Method: competitive inhibition

Note: This product is for in vitro research use only