

CHAC1 Polyclonal Antibody

Catalog No: #49244

Package Size: #49244-1 50ul #49244-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CHAC1 Polyclonal Antibody
Isotype	IgG
Purification	CHAC1 Polyclonal Antibody detects endogenous levels of protein.
Applications	WB;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	CHAC1 Polyclonal Antibody detects endogenous levels of protein.
Immunogen Description	Synthesized peptide derived from part region of human protein
Target Name	CHAC1 BOTCH
Uniprot	Q9BUX1
GeneID	79094
Calculated MW	24kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000 ELISA 1:5000-20000

Images



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only