

Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit



Catalog No: #EK20729

Orders: order@signalwayantibody.com

Package Size: #EK20729-1 48T #EK20729-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Alpha-Melanocyte Stimulating Hormone (aMSH) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:12.5-800pg/mL

Sensitivity:5.7pg/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only