

C1q-C Polyclonal Antibody

Catalog No: #49246

Package Size: #49246-1 50ul #49246-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	C1q-C Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	C1q-C Polyclonal Antibody detects endogenous levels of C1q-C protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human C1QC.
Target Name	C1QC
Other Names	C1QC;C1QG;Complement C1q subcomponent subunit C
Uniprot	P02747
GeneID	714
Calculated MW	30kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

Images



Western blot analysis of lysates from rat lung, using C1QC Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human lung carcinoma tissue, using C1QC Antibody. The picture on the right is blocked with the synthesized peptide.

Note: This product is for in vitro research use only