

GPX4 (Phospho-Ser104) Antibody

Catalog No: #SAB712P



Package Size: #SAB712P-1 50ul #SAB712P-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	GPX4 (Phospho-Ser104) Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antibodies were produced by immunizing rabbits with synthetic peptide and KLH conjugates. Antibodies were purified by affinity-chromatography using epitope-specific peptide.
Applications	Custom antibody
Species Reactivity	Hu
Immunogen Description	Immunogen Type:Peptide-KLH Immunogen Description:Peptide sequence around phosphorylation site of serine104 (C-AGYNVKFDMF-S(p)-KI) derived from Human GPX4.
Target Name	GPX4
Other Names	Phospholipid hydroperoxide glutathione peroxidase GPX4,PHGPx,Glutathione peroxidase 4 ,GPX4
Uniprot	P36969
GeneID	2879
Calculated MW	20/22kDa
Concentration	1.0mg/ml
Formulation	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg ²⁺ and Ca ²⁺), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C for long term preservation (recommended). Store at 4°C for short term use.

Note: This product is for in vitro research use only