

MFAP4 rabbit pAb

Catalog No: #49250



Package Size: #49250-1 50ul #49250-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MFAP4 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of MFAP4 at Human, Mouse
Immunogen Description	Synthesized peptide derived from human MFAP4
Target Name	MFAP4
Uniprot	P55083
GeneID	4239
Calculated MW	28kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Note: This product is for in vitro research use only