

F105B rabbit pAb

Catalog No: #49331

Package Size: #49331-1 50ul #49331-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	F105B rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of F105B at Human/Mouse
Immunogen Description	Synthesized peptide derived from human F105B
Target Name	FAM105B
Other Names	Protein FAM105B
Uniprot	Q96BN8
GeneID	90268
Calculated MW	38kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000

Images



Note: This product is for in vitro research use only