

Recombinant Human Polyadenylate-binding protein 2 (PABPN1)



Catalog No: #AP89684

Orders: order@signalwayantibody.com

Package Size: #AP89684-1 20ug #AP89684-2 100ug #AP89684-3 1mg

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human Polyadenylate-binding protein 2 (PABPN1)
Brief Description	Recombinant Protein
Other Names	Nuclear poly(A)-binding protein 1; OPMD; PAB2; PABII; PABP 2; pABP-2; PABP2; PABP2_HUMAN; PABPII; Pabpn1; poly(A) binding protein nuclear 1; Poly(A)-binding protein 2; Poly(A)-binding protein II; PolyA binding protein II; Polyadenylate-binding nuclear protein 1; Polyadenylate-binding protein 2
Uniprot	Q86U42
Target Species	Homo sapiens (Human)
Tag Info	Tag type will be determined during the manufacturing process.
Target Sequence	Protein Length: Full Length of Mature Protein Expression Region: 2-306 AAAAAAAAA AGAAGGRGSG PGRRRHLVPG AGGEAGEGAP GGAGDYGNGL ESEELEPEEL LLEPEPEPEP EEEPPRPRAP PGAPGPGPGS GAPGSQEEEE EPGLVEGDPG DGAIEDPELE AIKARVREME EEAELKELQ NEVEKQMNMS PPPGNAGPVI MSIEEKMEAD ARSIYVGNVD YGATAEELEA HFHGCOSVNR VTILCDKFSG HPGKFAYIEF SDKESVRTSL ALDESLFRGR QIKVIPKRTN RPGISTTDRG FPRARYRART TNYNSSRSRF YSGFNSRPRG RVYRGRARAT SWYSPY
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.

Application Details

>85% (SDS-PAGE)

Note: This product is for in vitro research use only