

MAP7D2 Rabbit pAb

Catalog No: #49348

Package Size: #49348-1 50ul #49348-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	MAP7D2 Rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	affinity purified by Protein A
Applications	WB, IHC-P, IHC-F, ICC, IF
Species Reactivity	Human;Mouse;Rat
Immunogen Description	KLH conjugated synthetic peptide derived from human MAP7D2: 631-732/732
Conjugates	Unconjugated
Target Name	MAP7D2
Other Names	FLJ14503; MAP7 domain containing 2; MAP7 domain containing 2 protein variant 1; MAP7 domain containing 2 protein variant 2; MAP7 domain-containing protein 2; RP11-393H10.2
Uniprot	Q96T17
GeneID	256714
Calculated MW	82 kDa
Concentration	1 mg/ml
Formulation	0.01M TBS (pH7.4) with 1% BSA, 0.02% Proclin300 and 50% Glycerol.
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.

Application Details

WB 1:500-2000 IHC-P 1:100-500 IHC-F 1:100-500 ICC 1:100-500 IF 1:100-500

Note: This product is for in vitro research use only