

# Porcine Immunoglobulin M (IgM) ELISA Kit

Catalog No: #EK20733



Package Size: #EK20733-1 48T #EK20733-2 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Porcine Immunoglobulin M (IgM) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Porcine

## Application Details

Detect Range: 0.312-20 µg/mL

Sensitivity: 0.103 µg/mL

Sample Type: serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 µL

Assay Time: 1-4.5h

Detection wavelength: 450 nm

## Product Description

Detection Method: sandwich

Note: This product is for in vitro research use only