

CD38 (ABT-CD38) mouse mAb (Ready to Use)

Catalog No: #58100

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	CD38 (ABT-CD38) mouse mAb (Ready to Use)
Brief Description	Mouse Monoclonal
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Human
Specificity	The antibody can specifically recognize human CD38 protein.
Immunogen Description	Synthesized peptide derived from human CD38
Target Name	CD38
Uniprot	P28907
GeneID	952
Formulation	The prediluted ready-to-use antibody is diluted in phosphate buffer saline containing stabilizing protein and 0.05% Proclin 300
Storage	2°C to 8°C/1 year

Application Details

Ready to use for IHC

Images



Human tonsil tissue was stained with Anti-CD38 (ABT-CD38) Antibody

Note: This product is for in vitro research use only