

S100 (ABT-S100B) mouse mAb (Ready to Use)

Catalog No: #58101

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	S100 (ABT-S100B) mouse mAb (Ready to Use)
Brief Description	Mouse Monoclonal
Host Species	Mouse
Clonality	Monoclonal
Isotype	Mouse IgG1, Kappa
Purification	The antibody was affinity-purified from ascites by affinity-chromatography using specific immunogen.
Applications	IHC
Species Reactivity	Human
Specificity	This antibody detects endogenous levels of S100 protein.
Immunogen Description	Synthesized peptide derived from human S100
Target Name	S100B
Uniprot	P04271
GeneID	6285
Formulation	The prediluted ready-to-use antibody is diluted in phosphate buffer saline containing stabilizing protein and 0.05% Proclin 300
Storage	2°C to 8°C/1 year

Application Details

Ready to use for IHC

Note: This product is for in vitro research use only