

## TRI58 rabbit pAb

Catalog No: #49640

Package Size: #49640-1 50ul #49640-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TRI58 rabbit pAb
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	This antibody detects endogenous levels of TRI58 at Human/Mouse
Immunogen Description	Synthesized peptide derived from human TRI58
Target Name	TRIM58
Uniprot	Q8NG06
GeneID	25893
Calculated MW	53kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1 $\mu$ g 500-2000

## Images



Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night

Note: This product is for in vitro research use only