

## IgA Polyclonal Antibody

Catalog No: #49714

Package Size: #49714-1 50ul #49714-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	IgA Polyclonal Antibody
Brief Description	rabbit polyclonal antibody
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Rat;Mouse
Specificity	IgA Polyclonal Antibody detects endogenous levels of IgA protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human IgA.
Target Name	IGHA1/IGHA2
Other Names	IGHA1;Ig alpha-1 chain C region;IGHA2;Ig alpha-2 chain C region
Uniprot	P01876/P01877
GeneID	3493/3494
Calculated MW	38kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000.. IF 1:50-200

## Images



Western Blot analysis of various cells using IgA Polyclonal Antibody



Immunohistochemistry analysis of IgA antibody in paraffin-embedded human lung carcinoma tissue.

---

Note: This product is for in vitro research use only