

CLK2 Polyclonal Antibody

Catalog No: #49996

Package Size: #49996-1 50ul #49996-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	CLK2 Polyclonal Antibody
Isotype	IgG
Applications	WB;IHC;IF;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	CLK2 Polyclonal Antibody detects endogenous levels of CLK2 protein.
Immunogen Description	The antiserum was produced against synthesized peptide derived from human CLK2.
Target Name	CLK2
Other Names	CLK2;Dual specificity protein kinase CLK2;CDC-like kinase 2
Uniprot	P49760
GeneID	1196
Calculated MW	60kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:20000.. IF 1:50-200

Images



Western blot analysis of lysates from COLO205 and A549 cells, using CLK2 Antibody. The lane on the right is blocked with the synthesized peptide.

Immunohistochemical analysis of paraffin-embedded Human brain.



Note: This product is for in vitro research use only