

Pacific Blue anti-mouse CD45

Catalog No: #26846



Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description

Product Name	Pacific Blue anti-mouse CD45
Isotype	Rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	Pacific Blue
Other Names	T200, Ly-5, LCA
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionoO PH7.2-7.4
Storage	Store at 4°C in dark

Application Details

FC 5μl/time

Note: This product is for in vitro research use only