

FITC anti-human/mouse CD11b

Catalog No: #26847



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	FITC anti-human/mouse CD11b
Isotype	rat IgG2b
Applications	FC
Species Reactivity	human/mouse
Conjugates	FITC
Other Names	α M integrin, Mac-1, Mo1, CR3, Ly-40, C3biR, ITGAM
Formulation	0.2%BSA ϵ 0.1%Proclin 300 Phosphate buffered solution \circ PH7.2-7.4
Storage	Store at 4°C in dark

Application Details

FC 5 μ l/time

Note: This product is for in vitro research use only