

# PE anti-mouse F4/80

Catalog No: #26848



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	PE anti-mouse F4/80
Isotype	rat IgG2b
Applications	FC
Species Reactivity	Mouse
Conjugates	PE
Other Names	EMR1, Ly71
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionO PH7.2-7.4
Storage	Store at 4°C in dark

## Application Details

FC 5&mu;l/time

Note: This product is for in vitro research use only