

Alexa Fluor 700 anti-mouse CD86

Catalog No: #26849



Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Alexa Fluor 700 anti-mouse CD86
Isotype	Rat IgG2a, k
Applications	FC
Species Reactivity	Mouse
Conjugates	Alexa Fluor 700
Other Names	B7-2, B70, Ly-58
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionO PH7.2-7.4
Storage	Store at 4°C in dark

Application Details

FC 5μl/time

Note: This product is for in vitro research use only