

# APC anti-mouse CD206

Catalog No: #26850



Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	APC anti-mouse CD206
Applications	FC
Species Reactivity	Mouse
Conjugates	APC
Formulation	0.2%BSAε 0.1%Proclin 300 Phosphate buffered solutionO PH7.2-7.4
Storage	Store at 4°C in dark

## Application Details

FC 5&mu;l/time

Note: This product is for in vitro research use only