

Human Methylcytosine dioxygenase TET3 (TET3) ELISA Kit

Catalog No: #EK20746

Orders: order@signalwayantibody.com

Package Size: #EK20746-1 48T #EK20746-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Methylcytosine dioxygenase TET3 (TET3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.781-50 ng/mL

Sensitivity:0.32ng/mL

Sample Type:serum, plasma, tissue homogenates, cell lysates, cell culture supernates or other biological fluids.

Sample Volume: 50 μ L

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only