

# Human SPDL1/ Protein Spindly ELISA Kit

Catalog No: #EK21137



Package Size: #EK21137-1 96T

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Human SPDL1/ Protein Spindly ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human
Other Names	Arsenite-related gene 1 protein, Coiled-coil domain-containing protein 99, Rhabdomyosarcoma antigen MU-RMS-40.4A, Spindle apparatus coiled-coil domain-containing protein 1, hSpindly, CCDC99
Uniprot	Q96EA4

## Application Details

Detect Range:0.312-20ng/mL

Sensitivity:0.083ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 &mu;L

Assay Time:1-4.5h

Detection wavelength:450 nm

## Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only