

Human Coiled-coil domain-containing protein 134 (CCDC134) ELISA Kit



Catalog No: #EK20749

Orders: order@signalwayantibody.com

Package Size: #EK20749-1 48T #EK20749-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Coiled-coil domain-containing protein 134 (CCDC134) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:0.156-10ng/mL
Sensitivity:0.067ng/mL
Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids
Sample Volume: 50 μL
Assay Time:1-4.5h
Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only