

Anti-Caspase 1 p20 Antibody

Catalog No: #62002

Package Size: #62002-1 50ul #62002-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

Product Name	Anti-Caspase 1 p20 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was purified by immunogen affinity chromatography.
Applications	WB
Species Reactivity	H, M, R
Immunogen Description	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human Caspase 1 p20. The exact sequence is proprietary.
Target Name	CASP1
Other Names	IL1BC; IL1BCE; Caspase-1; CASP-1; Interleukin-1 beta convertase; IL-1BC; Interleukin-1 beta-converting enzyme; ICE; IL-1 beta-converting enzyme; p45
Uniprot	P29466
GeneID	834
Formulation	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.

Application Details

WB (1/500 - 1/1000)

Images



Western blot analysis of Caspase 1 p20 expression in HEK293T UV-treated (A), Raw264.7 LPS-treated (B), rat heart (C) whole cell lysates.

Note: This product is for in vitro research use only