

Human Complement C1q Tumor Necrosis Factor Related protein 7 (C1QTNF7) ELISA Kit



Catalog No: #EK20750

Orders: order@signalwayantibody.com

Package Size: #EK20750-1 48T #EK20750-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Human Complement C1q Tumor Necrosis Factor Related protein 7 (C1QTNF7) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:28pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only