

Mouse 27-hydroxycholesterol (17-OHC) ELISA Kit

Catalog No: #EK20758



Package Size: #EK20758-1 48T #EK20758-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse 27-hydroxycholesterol (17-OHC) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:78.1-5000pg/mL

Sensitivity:27pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:competitive inhibition

Note: This product is for in vitro research use only