

Mouse Alkaline ceramidase 3 (ACER3) ELISA Kit

Catalog No: #EK20770



Package Size: #EK20770-1 48T #EK20770-2 96T

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Alkaline ceramidase 3 (ACER3) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:31.2-2000pg/mL

Sensitivity:12.5pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only