

Mouse Phosphatidylcholine:ceramide cholinephosphotransferase 1 (SGMS1) ELISA Kit

Catalog No: #EK20780

Orders: order@signalwayantibody.com

Package Size: #EK20780-1 48T #EK20780-2 96T

Support: tech@signalwayantibody.com

Description

Product Name	Mouse Phosphatidylcholine:ceramide cholinephosphotransferase 1 (SGMS1) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Mouse

Application Details

Detect Range:0.156-10ng/mL

Sensitivity:0.067ng/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only