

Human Interferon-inducible double stranded RNA-dependent protein kinase activator A (PRKRA) ELISA Kit

Catalog No: #EK20766

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Package Size: #EK20766-1 48T #EK20766-2 96T

Description

Product Name	Human Interferon-inducible double stranded RNA-dependent protein kinase activator A (PRKRA) ELISA Kit
Brief Description	ELISA Kit
Applications	ELISA
Species Reactivity	Human

Application Details

Detect Range:125-8000pg/ml

Sensitivity:43pg/mL

Sample Type:Serum, plasma, tissue homogenates, cell culture supernates and other biological fluids

Sample Volume: 50 μL

Assay Time:1-4.5h

Detection wavelength:450 nm

Product Description

Detection Method:sandwich

Note: This product is for in vitro research use only