

C/EBP  $\beta$  Monoclonal Antibody

Catalog No: #62014

Package Size: #62014-1 50ul #62014-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


## Description

Product Name	C/EBP $\beta$ Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	Affinity purification
Applications	WB;FCM;ELISA
Species Reactivity	Human
Specificity	C/EBP $\beta$ Monoclonal Antibody detects endogenous levels of C/EBP $\beta$ protein.
Immunogen Description	Purified recombinant fragment of human C/EBP $\beta$ expressed in E. Coli.
Target Name	CEBPB
Other Names	CEBPB;LAP;TCF5;PP9092;CCAAT/enhancer-binding protein beta;C/EBP beta;Liver activator protein;Nuclear factor NF-IL6;Transcription factor 5;TCF-5
Uniprot	P17676
GeneID	1051
Calculated MW	36kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

WB 1:500 - 1:2000. Flow cytometry: 1:200 - 1:400. ELISA: 1:10000.

## Images



Western Blot analysis using C/EBP  $\beta$  Monoclonal Antibody against human C/EBP  $\beta$  (AA: 161-338) recombinant protein.

Note: This product is for in vitro research use only