

## Biotin anti-human IL-12p40

Catalog No: #FC0147



Package Size: #FC0147-1 25ug #FC0147-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Biotin anti-human IL-12p40
Brief Description	human monoclonal antibody
Host Species	human
Clonality	Monoclonal
Clone No.	s-ABT-874
Isotype	human IgG1, $\lambda$
Applications	FC
Species Reactivity	human
Conjugates	Biotin
Target Name	IL-12p40
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5  $\mu$ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only