

PE-Cy5 anti-human/pig CD62E/CD62P

Catalog No: #FC0181



Package Size: #FC0181-1 25Tests #FC0181-2 100Tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


Product Name	PE-Cy5 anti-human/pig CD62E/CD62P
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-1.2B6
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human/pig
Conjugates	PE-Cy5
Target Name	CD62E/CD62P
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human Activated platelets stained with PE-Cy5 anti-human CD62E/CD62P (colored histogram) or PE-Cy5 mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only