

## Alexa Fluor 488 anti-human CD107b

Catalog No: #FC0360

Package Size: #FC0360-1 25Tests #FC0360-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

## Description


Product Name	Alexa Fluor 488 anti-human CD107b
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-H4B4
Isotype	mouse IgG1,k
Applications	FC
Species Reactivity	human
Conjugates	Alexa Fluor 488
Target Name	CD107b
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Camptothecin-stimulated cell line Jurkat was fixed permeabilized and intracellular stained with Alexa Fluor 488 anti-human CD107b (colored histogram) or Alexa Fluor 488 mouse IgG1 κ (gray histogram)

Note: This product is for in vitro research use only