

Alexa Fluor 700 anti-human/pig CD62E/CD62P

Catalog No: #FC0370



Package Size: #FC0370-1 25Tests #FC0370-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

Description


Product Name	Alexa Fluor 700 anti-human/pig CD62E/CD62P
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-1.2B6
Isotype	mouse IgG1
Applications	FC
Species Reactivity	human/pig
Conjugates	Alexa Fluor 700
Target Name	CD62E/CD62P
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 10 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human Activated platelets stained with Alexa Fluor 700 anti-human CD62E/CD62P (colored histogram) or Alexa Fluor 700 mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only