

# APC-Cy7 anti-human HLA-DR

Catalog No: #FC0461



Package Size: #FC0461-1 25Tests #FC0461-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)  
Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	APC-Cy7 anti-human HLA-DR
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-L243
Isotype	mouse IgG2a, κ
Applications	FC
Species Reactivity	human
Conjugates	APC-Cy7
Target Name	HLA-DR
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with APC-Cy7 anti-human HLA-DR (colored histogram) or APC-Cy7 mouse IgG2a κ (gray histogram)

Note: This product is for in vitro research use only