

PerCP-Cy5.5 anti-human CD38

Catalog No: #FC0691

Package Size: #FC0691-1 25Tests #FC0691-2 100Tests

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com

Description


| | |
|--------------------|---|
| Product Name | PerCP-Cy5.5 anti-human CD38 |
| Brief Description | mouse monoclonal antibody |
| Host Species | mouse |
| Clonality | Monoclonal |
| Clone No. | s-HIT2 |
| Isotype | mouse IgG1 |
| Applications | FC |
| Species Reactivity | human |
| Conjugates | PerCP-Cy5.5 |
| Target Name | CD38 |
| Formulation | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL. |

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PerCP-Cy5.5 anti-human CD38 (colored histogram) or PerCP-Cy5.5 mouse IgG1 (gray histogram)

Note: This product is for in vitro research use only