

# PerCP-Cy5.5 anti-human CD45RA

Catalog No: #FC0693

Package Size: #FC0693-1 25Tests #FC0693-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


Product Name	PerCP-Cy5.5 anti-human CD45RA
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-H1100
Isotype	mouse IgG2b
Applications	FC
Species Reactivity	human
Conjugates	PerCP-Cy5.5
Target Name	CD45RA
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Human peripheral blood lymphocytes stained with PerCP-Cy5.5 anti-human CD45RA (colored histogram) or PerCP-Cy5.5 mouse IgG2b (gray histogram)

Note: This product is for in vitro research use only