

PE-Cy7 anti-human CD42b

Catalog No: #FC0832



Package Size: #FC0832-1 25Tests #FC0832-2 100Tests

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	PE-Cy7 anti-human CD42b
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-HIP1
Isotype	mouse IgG1 α 0 Ω 1/2 Ω 1/2k
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy7
Target Name	CD42b
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 μ l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only