

PE-Cy7 anti-human CD57

Catalog No: #FC0845

Package Size: #FC0845-1 25Tests #FC0845-2 100Tests

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description


Product Name	PE-Cy7 anti-human CD57
Brief Description	mouse monoclonal antibody
Host Species	mouse
Clonality	Monoclonal
Clone No.	s-HNK-1
Isotype	mouse IgM
Applications	FC
Species Reactivity	human
Conjugates	PE-Cy7
Target Name	CD57
Formulation	0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4
Storage	Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.

Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10⁶ cells can be tested per experiment. Refer to the detailed scenario as you test.

Images



Human peripheral blood lymphocytes stained with PE-Cy7 anti-human CD57 (colored histogram) or PE-Cy7 mouse IgM (gray histogram)

Note: This product is for in vitro research use only