

# Biotin anti-human /mouse CD45R/B220

Catalog No: #FC0927



Package Size: #FC0927-1 25ug #FC0927-2 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                    |   |
|--------------------|---|
| Product Name       | Biotin anti-human /mouse CD45R/B220                               |
| Clonality          | Monoclonal  |
| Clone No.          | s-RA3-6B2   |
| Isotype            | rat IgG2a   |
| Applications       | FC  |
| Species Reactivity | human /mouse  |
| Conjugates         | Biotin  |
| Target Name        | CD45R/B220  |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 &mu;l/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

Note: This product is for in vitro research use only