

## FITC anti-mouse CD152

Catalog No: #FC1021

Package Size: #FC1021-1 25Tests #FC1021-2 100Tests

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description


|                    |   |
|--------------------|---|
| Product Name       | FITC anti-mouse CD152   |
| Brief Description  | Syrian Hamster Monoclonal Antibody                                |
| Host Species       | Syrian Hamster  |
| Clonality          | Monoclonal  |
| Clone No.          | s-9H10  |
| Isotype            | Syrian Hamster IgG  |
| Applications       | FC  |
| Species Reactivity | mouse   |
| Conjugates         | FITC  |
| Target Name        | CD152   |
| Formulation        | 0.2%BSA and 0.1%Proclin 300 Phosphate-buffered solution,PH7.2-7.4 |
| Storage            | Store at 4°C. DO NOT FREEZE. LIGHT SENSITIVE MATERIAL.            |

## Application Details

Test: 5 µl/Test

Notice: This reagent has been pre-diluted for use at recommended volume per test in flow cytometry analysis. Typically, 1x 10<sup>6</sup> cells can be tested per experiment. Refer to the detailed scenario as you test.

## Images



Note: This product is for in vitro research use only